Original CMT

California Mastitis Test

US Pat. Nos. 2998392, 2935384 and 3102418

Is mastitis widespread ?

Surprisingly enough, mastitis appears to be much more widespread than most realize. Dairy authorities have estimated that as many as one out of two cows has some mastitis in one or more quarters. Not all are serious cases, of course.

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What is CMT ?

CMT is an abbreviation for the original California Mastitis Test developed by Dr O. W. Schalm of the University of California's School of Veterinary Medicine for the detection of increased leukocyte cell count in fresh milk, resulting from mastitis. It's a simple test to make and interpret. CMT tells in seconds whether milk is mastitic and pinpoints the quarter delivering mastitic milk.

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How do you use CMT ?

The CMT Kit contains a bottle of CMT ready to use or concentrate (enough to make 1 gallon of test solution) and paddle. Complete instruction on diluting the concentrate are on the label. Detail on administering the CMT procedure are illustrated and described on this document.

Primary value of CMT is at the side of the cow where it aids in quick diagnosis of mastitis, as well as control, for it pinpoints the quarter secreting mastitic milk.

CMT also can be use to test bulk milk. When doing this, it's important to keep the milk cold and not to let it become too old. After 24 hours, aging and/or bacterial action tend to destroy factors responsible for positive reaction to the test.

Mastitic milk tends to jell when tested by the CMT procedure. The degree of jelling indicates the presence and severity of mastitis.

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Was the test positive ?

If CMT shows a positive mastitis reaction, as determined by consulting the chart in this document, call your veterinarian at once so that he can determine the extent and nature of the treatment necessary. Follow his instruction.

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After discarding the first stream of milk, draw the next milk into the shallow cups on the paddle, keeping the quarters separate.

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The ideal amount of milk is that which remains in the cup when the paddle is tilted to an almost vertical position, as shown.

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Form pools of milk in cups, as shown, by tilting paddle. Squirt test solution over milk. Avoid making bubbles. Proportion of solution to milk should be at least one to one.

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To mix solution and milk gently rotating paddle in horizontal plane, swirling the mixture. Positive reactions occur and can be graded during this rotary motion.

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Rhinehart Développement Corp.
Dairy Research Product, Inc.